

Meniscus

Dimensions

Description		W <i>cm</i>	A/AP <i>mm</i>	C/AP <i>mm</i>	WM/ML <i>mm</i>	WL/ML <i>mm</i>	TL <i>mm</i>	TM <i>mm</i>	DL <i>mm</i>	DM <i>mm</i>
530400	Meniscus with tibia plateau left	60.0-100.0	30.0-65.0	30.0-65.0	30.0-53.0	30.0-53.0	10.0-95.0	10.0-95.0	25.0-75.0	25.0-75.0
530410	Menisci with tibia plateau right	60.0-100.0	30.0-65.0	30.0-65.0	30.0-53.0	30.0-53.0	10.0-95.0	25.0-75.0	25.0-75.0	25.0-75.0
530402	Meniscus with hemi plateau lateral left	N/A	N/A	30.0-65.0	N/A	30.0-53.0	10.0-95.0	N/A	N/A	N/A
530412	Meniscus with hemi plateau lateral right	N/A	N/A	30.0-65.0	N/A	30.0-53.0	10.0-95.0	N/A	N/A	N/A
530401	Meniscus with hemi plateau medial left	N/A	30.0-65.0	N/A	30.0-53.0	N/A	N/A	10.0-95.0	N/A	N/A
530411	Meniscus with hemi plateau medial right	N/A	30.0-65.0	N/A	30.0-53.0	N/A	N/A	N/A	N/A	N/A
530501	Meniscus medial, isolated, left	N/A	30.0-65.0	N/A	30.0-53.0	N/A	N/A	N/A	N/A	N/A
530511	Meniscus medial, isolated, right	N/A	30.0-65.0	NA	30.0-53.0	N/A	N/A	N/A	N/A	N/A
530502	Meniscus lateral, isolated, left	N//A	N/A	30.3-65.0	N/A	30.0-53.0	N/A	N/A	N/A	N/A

530512	Meniscus lateral, isolated, right	N/A	N/A	30.0-65.0	N/A	30.0-53.0	N/A	N/A	N/A	N/A
--------	-----------------------------------	-----	-----	-----------	-----	-----------	-----	-----	-----	-----

This product is not in our standard product range. Please contact us for more information

Measurements

- W Maximum measure medial to lateral of the plateau including the bone
- DM Medial superior plateau measured anterior to posterior at the longest point
- DL Lateral superior plateau measured anterior to posterior at the longest point
- WM/ML Medial and lateral meniscus measured medial to lateral
- WL/ML Medial and lateral meniscus measured medial to lateral
- A/AP Medial and lateral meniscus measured anterior to posterior
- C/AP Medial and lateral meniscus measured anterior to posterior
- T Minimum thickness of the tibia not including soft tissue
- TM Minimum thickness of medial meniscus; do not release the attachment to the plateau for hemi
- TL Minimum thickness of the lateral meniscus-do not release the attachment to the plateau for hemi

